

PD1 antibody [RMP1-14] (Low endotoxin, azide free)

Cat. No. GTX05080

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA, Neutralizing /Inhibition
Reactivity	Mouse

Package
1 mg

PRODUCT

Summary Aggregation: <5% (Determined by SEC)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Neutralizing /Inhibition	Assay dependent
Not tested in other applications.	

Calculated MW 32 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Syrian Hamster BKH cells transfected with mouse PD-1 cDNA
Purification	Purified by Protein G From tissue culture supernatant
Purity	>95% by SDS-PAGE
Endotoxin	<0.001EU/μg (determined by the LAL assay)
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 09 / 02 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX05080 WB Image

WB analysis of various loading amounts of C-terminal His tagged mouse PD-1 proteins using GTX05080 PD1 antibody [RMP1-14] (Low endotoxin, azide free).

Lane 1 : 4μg

Lane 2 : 2μg

Dilution : 15μg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).