

CD86 antibody [C86/6500R]

Cat. No. GTX05083

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1-2 µg/ml

Note : Requires heating tissue sections in Tris-EDTA (pH 9.0)

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment (around aa 66-195) of human CD86 protein. The exact sequence is proprietary.
Purification	Purified by Protein A/G
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05083 IHC-P Image

IHC-P analysis of human lymph node tissue using GTX05083 CD86 antibody [C86/6500R].
 Antigen retrieval : Heat mediated antigen retrieval performed with Tris/EDTA buffer (pH 9.0)
 Dilution : 2 µg/ml


GTX05083 IHC-P Image

IHC-P analysis of human lymph node tissue using GTX05083 CD86 antibody [C86/6500R].
 Lower left panel: Negative control
 Antigen retrieval : Heat mediated antigen retrieval performed with Tris/EDTA buffer (pH 9.0)


GTX05083 IHC-P Image

IHC-P analysis of mouse tonsil tissue using GTX05083 CD86 antibody [C86/6500R].



For full product information, images and publications, please visit our [website](https://www.genetex.com).