

UGCG antibody

Cat. No. GTX05087

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1:5000
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 45 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	TBS, 1% BSA, 50% Glycerol
Preservative	0.02% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human UGCG.
Purification	Purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX05087 IHC-P Image

IHC-P analysis of rat small intestine tissue using GTX05087 UGCG antibody.

Antigen retrieval : Heat mediated antigen retrieval performed with sodium citrate buffer (pH 6.0)


Dilution : 1:200



GTX05087 WB Image

WB analysis of mouse small intestine tissue lysate using GTX05087 UGCG antibody.

Dilution : 1:1000



GTX05087 WB Image

WB analysis of mouse lymph node tissue lysate using GTX05087 UGCG antibody.

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).