

GNG2 antibody

Cat. No. GTX05089

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
200 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:8000
IHC-P	1:25

Not tested in other applications.

Calculated MW 8 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 19-52 amino acids from of human GNG2.
Purification	Protein A purified followed by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx05089 IHC-P Image

IHC-P analysis of human brain tissue using GTx05089 GNG2 antibody.

Antigen retrieval : Heat mediated antigen retrieval performed with citrate buffer (pH 6.0)

Dilution : 1:25



GTx05089 WB Image

WB analysis of various samples using GTx05089 GNG2 antibody.

Lane 1 : Human brain tissue lysate

Lane 2 : Mouse brain tissue lysate

Dilution : 1:8000

Loading : 20µg



For full product information, images and publications, please visit our [website](https://www.genetex.com).