

TAF15 antibody [27T75]

Cat. No. GTX05091

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ICC/IF	1:50-1:200
IHC-P	1:50-1:200
FCM	1:50

Not tested in other applications.

Calculated MW 62 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.4-0.5 mg/ml BSA, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthesized peptide derived from human TAF15.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05091 IHC-P Image

IHC-P analysis of human colorectal adenocarcinoma tissue using GTX05091 TAF15 antibody [27T75].

Antigen retrieval : Heat mediated antigen retrieval performed with EDTA buffer (pH 8.0)

Dilution : 1:50


GTX05091 WB Image

WB analysis of various samples using GTX05091 TAF15 antibody [27T75].

Dilution : 1:500


Loading : 30µg


GTX05091 IHC-P Image

IHC-P analysis of human cervix squamous cell carcinoma tissue using GTX05091 TAF15 antibody [27T75].

Antigen retrieval : Heat mediated antigen retrieval performed with EDTA buffer (pH 8.0)

Dilution : 1:50


GTX05091 IHC-P Image

IHC-P analysis of human prostate cancer tissue using GTX05091 TAF15 antibody [27T75].

Antigen retrieval : Heat mediated antigen retrieval performed with EDTA buffer (pH 8.0)

Dilution : 1:50



For full product information, images and publications, please visit our [website](https://www.genetex.com).