

PLA2R1 antibody

Cat. No. GTX05094

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	FCM, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	10µg/ml
ELISA	Assay dependent

Note : Permeabilization : 0.5% Triton
Antibody detection limit dilution : 1:1000

Not tested in other applications.

Calculated MW 169 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide derived from the internal region of human PLA2R1 protein (C-EKKKDTPKQHGT).
Purification	Purified by ammonium sulphate precipitation followed by antigen-affinity chromatography. From serum
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx05094 FCM Image

FCM analysis of PFA-fixed HeLa cells using GTx05094 PLA2R1 antibody.

Blue line : Primary antibody

Black line : Isotype control

Permeabilization : 0.5% Triton



For full product information, images and publications, please visit our [website](#).