

Haspin antibody

Cat. No. GTX05096

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, FCM, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.25-0.5 µg/ml
ICC/IF	5 µg/ml
FCM	1-3 µg/1x10 ⁶ fixed cells
ELISA	0.1-0.5 µg/ml

Note : Permeabilization step is needed.

Not tested in other applications.

Calculated MW 88 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E.coli-derived human Haspin recombinant protein (Position: R86-K798), which shares 65% amino acid (aa) sequence identity with mouse Haspin.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx05096 ICC/IF Image

ICC/IF analysis of PFA-fixed MG63 cells using GTx05096 Haspin antibody.

Antigen retrieval : Enzyme antigen retrieval reagent

Dilution : 5 µg/mL




GTx05096 WB Image

WB analysis of various samples using GTx05096 Haspin antibody.

Dilution : 0.5 µg/ml

Loading : 30µg



GTx05096 FCM Image

FCM analysis of 4% PFA-fixed K562 cells using GTx05096 Haspin antibody.

Blue line : Primary antibody

Green line : Isotype control

Red line : Blank

Dilution : 1 µg/1x10⁶ cells



For full product information, images and publications, please visit our [website](https://www.genetex.com).