

Human Calcium binding protein p22 protein, His tag

Cat. No. GTX05102-pro
Species Human

Package

50 µg

Applications
Observed MW (kDa) 27 kDa.

Properties
Form Lyophilized powder

Buffer Reconstitute with 10mM PBS to 0.1-1.0 mg/mL. Do not vortex. Lyophilized from PBS, 0.01% SKL, 5% Trehalose.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Region/Sequence N-terminal His tagged full length human calcium binding protein p22 protein (Gly2~His195, Uniprot #Q99653)

Expression System E. coli

Purity > 95% by SDS-PAGE


Endotoxin <1.0 EU/µg determined by LAL method

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX05102-pro Image

SDS-PAGE of GTX05102-pro Human Calcium binding protein p22 protein, His tag.


For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 10 / 10 Page 1 of 1