

Goat anti-Rabbit IgG antibody (TRITC)

Cat. No. GTX05104-25

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IHC-Fr, FCM, ELISA
Reactivity	Rabbit

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
ELISA	1:100-1:400

Not tested in other applications.

Product Note Preadsorbed with rabbit IgM. May react with immunoglobulins from other species.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit IgG
Purification	Purified by affinity chromatography on pooled rabbit IgG covalently linked to agarose From serum
Conjugation	Tetramethylrhodamine (TRITC) Wavelength


For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 12 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx05104-25 ELISA Image

ELISA analysis of rabbit IgG and IgM protein using serial diluted GTx05104-25 Goat anti-Rabbit IgG antibody (TRITC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).