

Protein A

Cat. No. GTX05108-pro

Applications WB, IP, ELISA, Purification

Package
5 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent
ELISA	Assay dependent
Purification	Assay dependent

Not tested in other applications.

Product Note

Binds the Fc portion of human IgG1, IgG2, and IgG4 as well as mouse IgG2a and IgG2b; also binds rabbit, porcine, guinea pig, and canine IgG.

Properties

Form Lyophilized powder

Buffer Reconstitute with PBS. Lyophilized from distilled water.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).