

CLK1 antibody

Cat. No. GTX05111

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ICC/IF	1:200-1:1000
IHC-P	1:100-1:300

Not tested in other applications.

Calculated MW 57 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.5% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthesized peptide derived from human CLK1. AA range:101-150
Purification	Purified by antigen-affinity chromatography From serum
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05111 WB Image

WB analysis of 293T whole cell lysate using GTX05111 CLK1 antibody.

Dilution : 1:500



GTX05111 IHC-P Image

IHC-P analysis of human lung cancer tissue using GTX05111 CLK1 antibody.

The picture on the right was treated with the blocking peptide.


Antigen retrieval : High pressure and heat mediated antigen retrieval performed with Tris-EDTA buffer (pH 8.0)

Dilution : 1:100


GTX05111 WB Image

WB analysis of 293T whole cell lysate using GTX05111 CLK1 antibody.

The lane on the right was treated with the blocking peptide.


GTX05111 WB Image

WB analysis of KB whole cell lysate using GTX05111 CLK1 antibody.

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).