

SUSD2 antibody [W5C5]

Cat. No. GTX05112

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
FCM	1-4 µg/ml

Note : Do not need a permeabilization step.

Not tested in other applications.

Calculated MW 90 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	15 mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	WERI-RB-1 retinoblastoma cells
Purification	Purified by Protein A
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05112 FCM Image

FCM analysis of surface staining MCF-7 cells using GTX05112 SUS D2 antibody [W5C5].

Red line : Primary antibody

Black line : Unstained cells

Dilution : 5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).