

SUSD2 antibody [W5C5] (APC)

Cat. No. GTX05112-07

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	1-4 µg/ml

Note : Do not need a permeabilization step.

Not tested in other applications.

Calculated MW 90 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	15 mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	WERI-RB-1 retinoblastoma cells
Purification	Purified by size-exclusion chromatography
Conjugation	Allophycocyanin (APC) Wavelength

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX05112-07 FCM Image

FCM analysis of surface staining HeLa cells using GTX05112-07 SUS2 antibody [W5C5] (APC).

Red line : Primary antibody

Black line : Isotype control

Dilution : 1.67 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).