

PI4 kinase beta antibody

Cat. No. GTX05121

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human, Mouse

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:50-1:100
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 91 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide 519-548 amino acids from the Central region of human PIK4CB.
Purification	Purified by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx05121 WB Image

WB analysis of HeLa whole cell lysate using GTx05121 PI4 kinase beta antibody.
Loading : 35µg



GTx05121 FCM Image

FCM analysis of HeLa cells using GTx05121 PI4 kinase beta antibody.
Green line : Primary antibody
Blue line : Negative control



GTx05121 IHC-P Image

IHC-P analysis of human skeletal muscle tissue using GTx05121 PI4 kinase beta antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).