

PML antibody

Cat. No. GTX05122

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 98 kDa. ([Note](#))

Product Note This antibody is specific for isoform 1, 9, 10 and 11 of PML protein.

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide directed towards the C terminal region of human PML: EGPSTLRVLDENLADPQAEDRPLVFFDLKIDNETQKISQLAAVNRESKFR
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX05122 IHC-P Image

IHC-P analysis of human spleen tissue using GTX05122 PML antibody.
Dilution : 5µg/ml




GTX05122 WB Image

WB analysis of HeLa whole cell lysate using GTX05122 PML antibody.
Dilution : 1µg/ml



GTX05122 IHC-P Image

IHC-P analysis of human bone marrow tissue using GTX05122 PML antibody.
Red : Primary antibody
Blue : DAPI
Dilution : 1:100



GTX05122 WB Image

WB analysis of human adult placenta tissue lysate using GTX05122 PML antibody.
Dilution : 1µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).