

Fibronectin antibody (ready-to-use)

Cat. No. GTX05128

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package
6 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent

Note : Recommendations for antigen retrieval : EDTA Buffer pH 8.0

Not tested in other applications.

Calculated MW 263 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Batch dependent
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Fibronectin isolated from a pool of normal human plasma.
Purification	Purified IgG From serum
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05128 IHC-P Image

IHC-P analysis of human kidney tissue using GTX05128 Fibronectin antibody (ready-to-use).



For full product information, images and publications, please visit our [website](#).