

Collagen IV antibody [COL-94] (ready-to-use)

Cat. No. GTX05132

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human

Package
6 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent

Note : Recommendations for antigen retrieval : HistoZyme (GTX03599)

Not tested in other applications.

Calculated MW 164 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Ascites fluid
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human collagen type IV.
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX05132 IHC-P Image

IHC-P analysis of human placenta tissue using GTX05132 Collagen IV antibody [COL-94] (ready-to-use).



For full product information, images and publications, please visit our [website](#).