

Hsp70 antibody [W27] (ready-to-use)

Cat. No. GTX05136

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	IHC-P
Reactivity	Human

Package
6 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent

Note : Recommendations for antigen retrieval : Citrate Buffer pH 6.0

Not tested in other applications.

Calculated MW	70 kDa. (Note)
Product Note	This antibody is specific to constitutive (HSP73) and inducible (HSP72) form of HSP70.

Properties

Form	Liquid
Buffer	Batch dependent
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	HSP70 from HeLa cells.
Purification	Purified IgG2a
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX05136 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX05136 Hsp70 antibody [W27] (ready-to-use).



For full product information, images and publications, please visit our [website](#).