

HBEGF antibody [HBEGF/9503]

Cat. No. GTX05155

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2-4 µg/ml
IHC-P	1-2 µg/ml

Note : Antigen retrieval : heat mediated tissue sections in 10mM Tris with 1mM EDTA (pH 9.0)

Not tested in other applications.

Calculated MW 23 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment (aa1-208) of human HBEGF protein
Purification	Purified by Protein A/G.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05155 WB Image

WB analysis of RAW 264.7 whole cell lysate using GTX05155 HBEGF antibody [HBEGF/9503].


GTX05155 IHC-P Image

IHC-P analysis of human skin tissue using GTX05155 HBEGF antibody [HBEGF/9503].

Lower left panel : Secondary antibody only



For full product information, images and publications, please visit our [website](https://www.genetex.com).