

## MAN1B1 antibody

Cat. No. GTX05162

|              |              |
|--------------|--------------|
| Host         | Rabbit       |
| Clonality    | Polyclonal   |
| Isotype      | IgG          |
| Applications | WB           |
| Reactivity   | Human, Mouse |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                | Recommended dilution |
|-----------------------------------|----------------------|
| WB                                | 1:500-1:2000         |
| Not tested in other applications. |                      |

Calculated MW 80 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 0.5% BSA, 50% Glycerol                                                                                                                                                                                                |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Synthesized peptide derived from MAN1B1 at AA range: 100-180                                                                                                                                                               |
| Purification  | Purified by antigen-affinity chromatography<br>From serum                                                                                                                                                                  |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx05162 WB Image

WB analysis of HeLa whole cell lysate using GTx05162 MAN1B1 antibody. The lane on the right was treated with blocking peptide.



### GTx05162 WB Image

WB analysis of 3T3 whole cell lysate using GTx05162 MAN1B1 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).