

Fas Ligand antibody

Cat. No. GTX05180

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1:5000
ICC/IF	1:50-1:200
IHC-P	1:200-1:400
FCM	1:20-1:100

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.01M TBS, 1% BSA, 50% Glycerol
Preservative	0.02% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human Fas Ligand
Purification	Purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTx05180 WB Image

WB analysis of Jurkat whole cell lysate using GTx05180 Fas Ligand antibody.

Dilution : 1:1000




GTx05180 IHC-P Image

IHC-P analysis of mouse placenta tissue using GTx05180 Fas Ligand antibody.

Antigen retrieval : Heat mediated antigen retrieval performed with sodium citrate buffer (pH 6.0).

Dilution : 1:200




GTx05180 WB Image

WB analysis of various samples using GTx05180 Fas Ligand antibody.

Lane 1 : Mouse bone tissue lysate

Lane 2 : Rat bone tissue lysate

Dilution : 1:1000



GTx05180 IHC-P Image

IHC-P analysis of human rectal carcinoma tissue using GTx05180 Fas Ligand antibody.

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).