

Succinyl Lysine antibody [ARC55837]

Cat. No. GTX05185

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, Dot
Reactivity	Species independent

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:5000
Dot	1:500-1:2000

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA, 50% Glycerol
Preservative	Proclin300 or sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence containing succinylated K.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05185 WB Image

WB analysis of mouse liver tissue lysate using GTX05185 Succinyl Lysine antibody [ARC55837].

Loading : 25µg


Dilution : 1:2000


GTX05185 Dot Image

Dot blot analysis of various peptides using GTX05185 Succinyl Lysine antibody [ARC55837].

Coating : 10 or 100 ng

Dilution : 1:1000


GTX05185 WB Image

WB analysis of HeLa cells with/without 100 nM Disodium succinate treatment using GTX05185 Succinyl Lysine antibody [ARC55837].

Loading : 30µg

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).