

Grp78 antibody

Cat. No. GTX05188

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Cryphonectria parasitica, Human, Mouse, Rat, Rabbit, Bovine, Dog, Hamster, Monkey, Fungi, Xenopus

Package
200 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:100
IHC-P	Assay dependent
IP	Assay dependent

Not tested in other applications.

Calculated MW 72 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rat GRP78 (Bip) synthetic peptide conjugated to KLH
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx05188 WB Image

WB analysis of various samples using GTx05188 Grp78 antibody.

Loading : 15µg

Dilution : 1:1000




GTx05188 IHC-P Image

IHC-P analysis of mouse brain tissue using GTx05188 Grp78 antibody.

Red : Primary antibody

Blue : DAPI

Dilution : 1:100



GTx05188 ICC/IF Image


ICC/IF analysis of 2% PFA fixed heat-shocked HeLa cells using GTx05188 Grp78 antibody.

Green : Primary antibody

Blue : DAPI

Dilution : 1:100

Magnification : 100x



GTx05188 WB Image

WB analysis of various samples using GTx05188 Grp78 antibody.

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).