

CPSF73 antibody [ARC2567]

Cat. No. GTX05201

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:4000
IP	0.5-4µg for 200-400µg extracts of whole cells

Not tested in other applications.

Calculated MW 77 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA, 50% Glycerol
Preservative	Proclin300 or sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthesized peptide (585-684 a.a.) derived from human CPSF3.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05201 WB Image

WB analysis of various samples using GTX05201 CPSF73 antibody [ARC2567].

Loading : 25 µg


Dilution : 1:1000


GTX05201 IP Image

IP analysis of HeLa whole cell lysate using GTX05201 CPSF73 antibody [ARC2567].

IP loading amount : 0.5 µg / 200µg extracts

Dilution : 1:1000


GTX05201 WB Image

WB analysis of rat spleen tissue lysate using GTX05201 CPSF73 antibody [ARC2567].

Loading : 25 µg

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).