

Fas Ligand antibody [MFL3]

Cat. No. GTX05207

Host	Armenian Hamster
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P, FCM, in vitro, in vivo, Neutralizing /Inhibition
Reactivity	Mouse

Package
1 mg

Applications

Calculated MW 31 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.2 µm filtered PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	BHK cells expressing B6 mouse FasL
Purification	Purified by Protein A. tissue culture supernatant
Purity	≥95% (Determined by SDS-PAGE)
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).