

Mycoplasma pneumoniae p116 protein, His tag

Cat. No. GTX05215-pro

Species	Mycoplasma pneumoniae
---------	-----------------------

Package	100 µg
---------	--------

Properties

Form	Liquid
------	--------

Buffer	PBS, 25mM arginine.
--------	---------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. For long-term storage, aliquot and store at -70°C or below. Do not vortex.
---------	--

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence	The Mycoplasma Pneumoniae P116 protein (protein accession number: CAA96035) fused to a His-Tag.
-----------------	---

Expression System	E.coli
-------------------	--------

Purification	Purified by GSH affinity chromatography
--------------	---

Purity	>95% by SDS-PAGE
--------	------------------

Conjugation	Unconjugated
-------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).