

Mycoplasma pneumoniae p30 protein, His tag

Cat. No. GTX05216-pro

Species Mycoplasma pneumoniae

Package
100 µg

Properties

Form Liquid

Buffer PBS, 25mM K₂CO₃

Preservative 0.025% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence Mycoplasma Pneumoniae P30 protein (protein accession number: P75330) is fused to a 6 amino acid His-tag at C-terminus.

Expression System E.coli

Purification Purified by chromatography technique.

Purity >95% by SDS-PAGE

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).