

KIF13A antibody

Cat. No. GTX05224

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-2000
ICC/IF	1:50-200

Not tested in other applications.

Calculated MW 202 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% glycerol.
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human KIF13A protein fragment(1300-1635 aa).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX05224 ICC/IF Image

ICC/IF analysis of L929 cells using GTX05224 KIF13A antibody.

Orange : Primary antibody

Blue : DAPI

Dilution : 1:100




GTX05224 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX05224 KIF13A antibody.

Orange : Primary antibody

Blue : DAPI

Dilution : 1:100



GTX05224 WB Image

WB analysis of LO2 whole cell lysate using GTX05224 KIF13A antibody.

Loading : 25 µg

Dilution : 1:1000



For full product information, images and publications, please visit our [website](#).