

TACSTD2 antibody [6H10B11]

Cat. No. GTX05229

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-2000
IHC-P	1:200-1000
FCM	1:200-400
ELISA	1:10000

Not tested in other applications.

Calculated MW 36 kDa. ([Note](#))

Product Note This antibody was raised against human TACSTD2 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human TACSTD2 (27-274 a.a.)
Purification	Purified IgG2b
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX05229 FCM Image

FCM analysis of HepG2 cells using GTX05229 TACSTD2 antibody [6H10B11].


Green : Primary antibody

Red : Negative control



GTX05229 IHC-P Image

IHC-P analysis of human cervical cancer tissue using GTX05229 TACSTD2 antibody [6H10B11].




GTX05229 FCM Image

FCM analysis of BEL-7402 cells using GTX05229 TACSTD2 antibody [6H10B11].

Green : Primary antibody

Red : Negative control



GTX05229 ELISA Image

ELISA analysis of antigen using serially diluted GTX05229 TACSTD2 antibody [6H10B11].

Black : Control antigen 100 ng

Purple : Antigen 10 ng

Blue : Antigen 50 ng

Red : Antigen 100 ng



For full product information, images and publications, please visit our [website](http://www.genetex.com).