

BACH1 antibody [ARC56290]

Cat. No. GTX05233

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ChIP assay
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:100-8000
ChIP assay	5µg for 10µg-15µg of Chromatin

Not tested in other applications.

Calculated MW 82 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.5% BSA, 50% Glycerol
Preservative	Proclin300 or sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant fragment of human BACH1 (410-736 a.a.)
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX05233 WB Image

Simple WB analysis of K562 whole cell lysate using GTX05233 BACH1 antibody [ARC56290].


Left : the target band (as indicated) in different concentrations

Right : chemiluminescence intensity against molecular weight along the capillary in different concentrations

Dilution : 1:100


GTX05233 ChIP assay Image


ChIP analysis of K562 whole cell lysate using GTX05233 BACH1 antibody [ARC56290].


GTX05233 WB Image

WB analysis of various samples using GTX05233 BACH1 antibody [ARC56290].

Loading : 25 μg

Dilution : 1:11000


GTX05233 WB Image

Simple WB analysis of mouse thymus tissue lysate using GTX05233 BACH1 antibody [ARC56290].

Left : the target band (as indicated) in different concentrations

Right : chemiluminescence intensity against molecular weight along the capillary in different concentrations

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).