

Streptavidin antibody [S3E11]

Cat. No. GTX10021

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2
Application	WB, FACS, ELISA, IHC
Reactivity	Streptomyces avidinii

Package
500 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FACS	Assay dependent
ELISA	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Product Note

No cross-reactivity with egg white avidin. Affinity (Kd) > 1 x 10⁸. Clone S8C12 reacts with the determinant nonoverlapping with clone S8E4, S10D4 and S3E11 epitopes. This antibody can bind to the streptavidin without interfering with the biotin-binding reaction.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Streptavidin
Purification	Protein A purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).