

FOXE1 antibody

Cat. No. GTX100275

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human FOXE1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx100275 IP Image

FOXE1 antibody immunoprecipitates FOXE1 protein in IP experiments.

IP Sample: HeLa whole cell lysate/extract

A. 40 µg HeLa whole cell lysate/extract


B. Control with 2 µg of preimmune rabbit IgG

C. Immunoprecipitation of FOXE1 protein by 2 µg of FOXE1 antibody (GTx100275)

15% SDS-PAGE

The immunoprecipitated FOXE1 protein was detected by FOXE1 antibody (GTx100275) diluted at 1:1000.

EasyBlot anti-rabbit IgG (GTx221666-01) was used as a secondary reagent.



GTx100275 WB Image

Sample (30 ug of whole cell lysate)

A: HeLaS3

10% SDS PAGE

GTx100275 diluted at 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).