

## D-Dimer antibody [DD3]

Cat. No. GTX10053

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2b
<b>Applications</b>	WB, ELISA, Immunoassay, Lateral Flow
<b>Reactivity</b>	Mouse

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Immunoassay	Assay dependent
Lateral Flow	Assay dependent

**Note : Immunoassays for the quantitative determination of D-dimer and high molecular weight fibrin degradation products.**

**Capture : GTX10053, Detection : GTX05175**

**Please notice that the detected antibody needs to be conjugated to Gold to function as the detection antibody when paired with the capture antibody.**

**Please contact us for custom Gold-conjugated antibody.**

Not tested in other applications.

**Product Note** This antibody recognize D-dimer and high molecular weight fibrin degradation products. It does not cross-react with fibrinogen.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	High molecular weight fibrin degradation products or synthetic peptides covering the cross-linked region of D-dimer gamma-chain.
<b>Purification</b>	Protein A purified



For full product information, images and publications, please visit our [website](#).

**Conjugation**

Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).