

Twist1 antibody [N1C3]

Cat. No. GTX100621

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:5000-1:20000
Not tested in other applications.	

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Twist1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX100621 WB Image


Twist1 antibody [N1C3] detects Twist1 protein by Western blot analysis.

A. 30 µg HeLa whole cell lysate/extract

B. 30 µg whole cell lysate/extract of Flag-human TWIST1-transfected HeLa cells

12 % SDS-PAGE

Twist1 antibody [N1C3] (GTX100621) dilution: 1:10000



GTX100621 WB Image

The data was published in the journal PLoS One in 2017. [PMID: 28683123](https://doi.org/10.1371/journal.pone.0186831)



For full product information, images and publications, please visit our [website](https://www.genetex.com).