

PB1 antibody [N1N2], N-term

Cat. No. GTX100781

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 193 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human PB1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx100781 IHC-P Image

PB1 antibody [N1N2], N-term detects PB1 protein at nucleus in mouse duodenum by immunohistochemical analysis.

Sample: Paraffin-embedded mouse duodenum.

PB1 antibody [N1N2], N-term (GTx100781) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTx100781 WB Image

Sample (30 µg of whole cell lysate)


A: A431 (GTx27909)

5% SDS PAGE

PB1 antibody

GTx100781 diluted at 1:500

The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.


GTx100781 ICC/IF Image

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using PB1(GTx100781) antibody at 1:500 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).