

DYNC1H1 antibody

Cat. No. GTX101452

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 532 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.89 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein corresponding to human DYNC1H1
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX101452 WB Image

DYNC1H1 antibody detects DYNC1H1 protein by western blot analysis. Whole cell extracts (30 μ g) was separated by 5% SDS-PAGE, and the membrane was blotted with DYNC1H1 antibody (GTX101452) at a dilution of 1:500.


GTX101452 WB Image

DYNC1H1 antibody detects DYNC1H1 protein by western blot analysis. Various whole cell extracts (30 μ g) were separated by 5 % SDS-PAGE, and blotted with DYNC1H1 antibody (GTX101452) diluted by 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).