

p62Dok antibody [N1C1]

Cat. No. GTX101610

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human p62Dok. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX101610 WB Image

Sample (30 ug of whole cell lysate)
A: Raji
10% SDS PAGE
GTX101610 diluted at 1:1000



GTX101610 IHC-P Image

Immunohistochemical analysis of paraffin-embedded H661 xenograft, using DOK1(GTX101610) antibody at 1:500 dilution.
Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX101610 WB Image

The data was published in the journal Int J Mol Sci in 2018. [PMID: 29518977](https://pubmed.ncbi.nlm.nih.gov/29518977/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).