

RASSF7 antibody

Cat. No. GTX101740

Host	Rabbit
Clonality	Polyclonal
Isotype	lgG
Applications	WB
Reactivity	Human

Package $100~\mu\text{l},\,25~\mu\text{l}$

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
Not tested in other applications.	

Calculated MW 40 kDa. (<u>Note</u>)

Properties	
Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human RASSF7. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our website.

Date 2025 / 12 / 13 Page 1 of 2



DATA IMAGES



GTX101740 WB Image

Sample (30 ug of whole cell lysate) A: Raji

B: K562

10% SDS PAGE

GTX101740 diluted at 1:1000



For full product information, images and publications, please visit our website.

Date 2025 / 12 / 13 Page 2 of 2