

RPL15 antibody

Cat. No. GTX101830

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 24 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.63 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human RPL15. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTx101830 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human normal colon, using RPL15(GTx101830) antibody at 1:500 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx101830 WB Image


Sample (30 ug of whole cell lysate)

A: H1299

12% SDS PAGE

RPL15 antibody

GTx101830 diluted at 1:1000




GTx101830 IHC-P Image

RPL15 antibody detects RPL15 protein at cytosol on human breast carcinoma by immunohistochemical analysis.

Sample: RPL15 antibody (GTx101830) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx101830 WB Image

The data was published in the journal Cell Death Dis in 2016. [PMID: 27253413](https://pubmed.ncbi.nlm.nih.gov/27253413/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).