

PSA antibody [PS1]

Cat. No. GTX10188

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Sandwich ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX10188/GTX10189, Detection : GTX10189/GTX10188

Not tested in other applications.

Product Note This antibody reacts with free PSA (epitope 1) . It does not cross react with kallikrein 2.

Properties

Form	Liquid
Buffer	10mM Tris, 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human PSA antigen purified from seminal plasma.
Purification	Purified by ion exchange chromatography
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).