

DNAJB6 antibody [N3C3]

Cat. No. GTX101947

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-Wm
Reactivity	Human, Zebrafish

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-Wm	Assay dependent

Not tested in other applications.

Calculated MW 36 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	0.1M Tris, 0.1M Glycine, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.58 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human DNAJB6. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX101947 IHC-Wm Image

Immunohistochemical analysis (whole mount) of zebrafish embryo, using DNAJB6 antibody [N3C3] (GTX101947) at 1:200 dilution.


GTX101947 WB Image

Whole cell extract (30 µg) was separated by 12% SDS-PAGE, and the membrane was blotted with DNAJB6 antibody [N3C3] (GTX101947) diluted at 1:500.


GTX101947 ICC/IF Image


DNAJB6 antibody [N3C3] detects DNAJB6 protein at nucleus by immunofluorescent analysis.

Sample: NT2D1 cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: DNAJB6 protein stained by DNAJB6 antibody [N3C3] (GTX101947) diluted at 1:500.

Blue: Hoechst 33342 staining.

Scale bar = 10 µm.


GTX101947 WB Image

Whole cell extract (30 µg) was separated by 12% SDS-PAGE, and the membrane was blotted with DNAJB6 antibody [N3C3] (GTX101947) diluted at 1:500.



For full product information, images and publications, please visit our [website](https://www.genetex.com).