

PAPP A antibody [18A10]

Cat. No. GTX10196

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 181 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human PAPPA antigen
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX10196 WB Image

Detection of human PAPP-A and proMBP subunits of htPAPP-A by monoclonal antibodies in Western blot. Lane 1: PAPP A antibody [18A10] (GTX10196); Lane 2: PAPP A antibody [3C8] (GTX10194); Lane 3: PAPP A antibody [5H9] (GTX10191); Lane 4: PAPP A antibody [11E4] (GTX14597); Lane 5: PAPP A antibody [10A5] (GTX10193); Lane 6: PAPP A antibody [3C8] (GTX10194); Lanes 1-4: reducing conditions; Lanes 5-6: non-reducing conditions, heterotetrameric complex was detected by PAPP A antibody [10A5] (GTX10193) and PAPP A antibody [3C8] (GTX10194)



For full product information, images and publications, please visit our [website](#).