

EIF3F antibody [N1C3]

Cat. No. GTX101967

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 38 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 20% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.39 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human EIF3F. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx101967 ICC/IF Image

Immunofluorescence analysis of methanol-fixed HeLa, using EIF3F(GTx101967) antibody at 1:200 dilution.


GTx101967 IHC-P Image

Immunohistochemical analysis of paraffin-embedded BT483 xenograft, using EIF3F(GTx101967) antibody at 1:100 dilution.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min



GTx101967 IHC-P Image

EIF3F antibody [N1C3] detects EIF3F protein at nucleus on mouse intestine by immunohistochemical analysis.

Sample: Paraffin-embedded mouse intestine.

EIF3F antibody [N1C3] (GTx101967) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTx101967 WB Image

Sample (30 ug of whole cell lysate)

A: A431

B: HepG2

10% SDS PAGE

GTx101967 diluted at 1:1000



For full product information, images and publications, please visit our [website](http://www.genetex.com).