

HAX1 antibody

Cat. No. GTX101992

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 32 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.43 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human HAX1. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX101992 IHC-P Image

HAX1 antibody detects HAX1 protein at cytoplasm in rat ovary by immunohistochemical analysis.

Sample: Paraffin-embedded rat ovary.


HAX1 antibody (GTX101992) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX101992 IHC-P Image

Immunohistochemical analysis of paraffin-embedded NCIN87 Xenograft , using HAX1(GTX101992) antibody at 1:100 dilution.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX101992 ICC/IF Image


HAX1 antibody detects HAX1 protein at cell membrane by immunofluorescent analysis.

Sample: HeLa cells were fixed in ice-cold MeOH for 5 min.

Green: HAX1 stained by HAX1 antibody (GTX101992) diluted at 1:500.

Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (GTX628802) diluted at 1:1000.

Blue: Fluoroshield with DAPI (GTX30920).


GTX101992 ICC/IF Image

Confocal immunofluorescence analysis (Olympus FV10i) of methanol-fixed A431, using HAX1(GTX101992) antibody (Green) at 1:500 dilution. Alpha-tubulin filaments were labeled with GTX11304 (Red) at 1:2000.



For full product information, images and publications, please visit our [website](https://www.genetex.com).