

GABARAPL2 antibody

Cat. No. GTX102006

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Wm, IHC
Reactivity	Human, Mouse, Rat, Zebrafish

References (1)

★★★★☆ Review (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000
IHC-Wm	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 14 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 40% Glycerol
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human GABARAPL2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx102006 WB Image

GABARAPL2 antibody detects GABARAPL2 protein by western blot analysis. Various whole cell extracts (30 µg) were separated by 15% SDS-PAGE, and the membrane was blotted with GABARAPL2 antibody (GTx102006) diluted at a dilution of 1:1500.



GTx102006 IHC Image

GABARAPL2 antibody detects Gabarapl2 protein on photoreceptors by immunohistochemical analysis. Sample: Agarose-embedded zebrafish embryo GABARAPL2 antibody (GTx102006) dilution: 1:200. Image provided with permission courtesy of Dr. T. Schilling at UC, Irvine.



GTx102006 IHC-Wm Image

GABARAPL2 antibody detects Gabarapl2 protein on zebrafish by whole mount immunohistochemical analysis.

Sample: 2 days-post-fertilization zebrafish embryo.

GABARAPL2 antibody (GTx102006) dilution: 1:50.



For full product information, images and publications, please visit our [website](https://www.genetex.com).