

TNF alpha antibody [1D9]

Cat. No. GTX10204

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC, Neutralizing /Inhibition
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC	Assay dependent
Neutralizing /Inhibition	Assay dependent

Not tested in other applications.

Product Note This antibody is effective in inhibition of the biological activity of recombinant and natural TNF- α .

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human recombinant tumor necrosis factor of alpha type.
Purification	Protein A purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).