

ADAMTSL2 antibody [C3], C-term

Cat. No. GTX102069

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (9)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 105 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human ADAMTSL2. The exact sequence is proprietary.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX102069 IHC-P Image

ADAMTSL2 antibody [C3], C-term detects ADAMTSL2 protein at cytosol(secreted) on SkHep1xenograft by immunohistochemical analysis.


Sample: Paraffin-embedded SkHep1xenograft.

ADAMTSL2 antibody [C3], C-term (GTX102069) dilution: 1:500.

Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min


GTX102069 WB Image

HepG2 whole cell extract and conditioned medium (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with ADAMTSL2 antibody [C3], C-term (GTX102069) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX102069 IHC-P Image

ADAMTSL2 antibody [C3], C-term detects ADAMTSL2 protein by immunohistochemical analysis.

Sample: Paraffin-embedded rat tissues.

ADAMTSL2 stained by ADAMTSL2 antibody [C3], C-term (GTX102069) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


GTX102069 WB Image

C2C12 whole cell and differentiated extracts (30 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with ADAMTSL2 antibody [C3], C-term (GTX102069) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.



For full product information, images and publications, please visit our [website](https://www.genetex.com).