

Cardiac Troponin T antibody [1F11]

Cat. No. GTX10214

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Application	Purification
Reactivity	Human, Mouse, Rat, Rabbit, Goat, Bovine, Cat, Dog, Pig, Fish

Package
200 µg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Purification	Assay dependent
Not tested in other applications.	
Calculated MW	36 kDa. (Note)
Product Note	No cross-reactivity with skeletal TnT.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Free human cTnT
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).