

Cardiac Troponin T antibody [7A9]

Cat. No. GTX10221

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, Immunoassay, Purification
Reactivity	Human, Rabbit, Goat, Bovine, Dog

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Immunoassay	Assay dependent
Purification	Assay dependent

Not tested in other applications.

Calculated MW 36 kDa. ([Note](#))

Product Note Does not cross-react with skeletal muscle troponin T, cTnI and TnC.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Full length protein (Human).
Purification	Protein A purified
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).