

Cardiac Troponin T antibody [7G7]

Cat. No. GTX10223

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Immunoassay, Purification, Sandwich ELISA
Reactivity	Human, Bovine

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Immunoassay	Assay dependent
Purification	Assay dependent
Sandwich ELISA	Assay dependent

Note : Best capture antibodies for BIAcore.**Capture : GTX10223, Detection : GTX28296/GTX28295/GTX10216**

Not tested in other applications.

Product Note No cross-reaction with skeletal TnT.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Free human cTnT
Purification	Protein A purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).